



# AIMP1 Antibody, FITC conjugated

|                            |                                                                                                                                                                                                                                                                                                                             |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA618638EC01HU                                                                                                                                                                                                                                                                                                          |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                                                                                                               |
| <b>Uniprot No.</b>         | Q12904                                                                                                                                                                                                                                                                                                                      |
| <b>Immunogen</b>           | Recombinant Human Aminoacyl tRNA synthase complex-interacting multifunctional protein 1 protein (25-180AA)                                                                                                                                                                                                                  |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                                                                                                                                                                      |
| <b>Species Reactivity</b>  | Human                                                                                                                                                                                                                                                                                                                       |
| <b>Form</b>                | Liquid                                                                                                                                                                                                                                                                                                                      |
| <b>Conjugate</b>           | FITC                                                                                                                                                                                                                                                                                                                        |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                                                                                                                                                                                                                                            |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                                                                                                                                                                                                                                    |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                                                                                                                         |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                                                                                  |
| <b>Alias</b>               | Aminoacyl tRNA synthase complex-interacting multifunctional protein 1 (Multisynthase complex auxiliary component p43) [Cleaved into: Endothelial monocyte-activating polypeptide 2 (EMAP-2) (Endothelial monocyte-activating polypeptide II) (EMAP-II) (Small inducible cytokine subfamily E member 1)], AIMP1, EMAP2 SCYE1 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                        |
| <b>Research Area</b>       | Cell Biology                                                                                                                                                                                                                                                                                                                |
| <b>Target Names</b>        | AIMP1                                                                                                                                                                                                                                                                                                                       |