



PLEC Antibody, FITC conjugated

Product Code	CSB-PA619965LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15149
Immunogen	Recombinant Human Plectin protein (1744-1929AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plectin (PCN) (PLTN) (Hemidesmosomal protein 1) (HD1) (Plectin-1), PLEC, PLEC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLEC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.