



ARFRP1 Antibody, HRP conjugated

Product Code	CSB-PA621696LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13795
Immunogen	Recombinant Human ADP-ribosylation factor-related protein 1 protein (1-201AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-ribosylation factor-related protein 1 (ARF-related protein 1) (ARP), ARFRP1, ARP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARFRP1