



LAPTM4A Antibody, Biotin conjugated

Product Code	CSB-PA621858LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15012
Immunogen	Recombinant Human Lysosomal-associated transmembrane protein 4A protein (181-233AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysosomal-associated transmembrane protein 4A (Golgi 4-transmembrane-spanning transporter MTP), LAPTM4A, KIAA0108 LAPTM4 MBNT MTRP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LAPTM4A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.