



PICALM Antibody, Biotin conjugated

Product Code	CSB-PA622660LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13492
Immunogen	Recombinant Human Phosphatidylinositol-binding clathrin assembly protein (517-624AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol-binding clathrin assembly protein (Clathrin assembly lymphoid myeloid leukemia protein), PICALM, CALM
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PICALM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.