



RPE65 Antibody, HRP conjugated

Product Code	CSB-PA624103LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16518
Immunogen	Recombinant Human Retinoid isomerohydrolase protein (348-484AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinoid isomerohydrolase (EC 3.1.1.64) (All-trans-retinyl-palmitate hydrolase) (Meso-zeaxanthin isomerase) (EC 5.3.3.-) (Retinal pigment epithelium-specific 65 kDa protein) (Retinol isomerase), RPE65
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	RPE65
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.