



PAG2 Antibody, HRP conjugated

Product Code	CSB-PA637747LB01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q28057
Immunogen	Recombinant Bovine Pregnancy-associated glycoprotein 2 protein (22-376AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pregnancy-associated glycoprotein 2 (PAG 2) (EC 3.4.23.-), PAG2
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	PAG2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.