



NADK2 Antibody, Biotin conjugated

Product Code	CSB-PA676825LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4G0N4
Immunogen	Recombinant Human NAD kinase 2, mitochondrial protein (16-156AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NAD kinase 2, mitochondrial (EC 2.7.1.23) (Mitochondrial NAD kinase) (NAD kinase domain-containing protein 1, mitochondrial), NADK2, C5orf33 MNADK NADKD1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NADK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.