



NEK5 Antibody, FITC conjugated

Product Code	CSB-PA740713LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P3R8
Immunogen	Recombinant Human Serine/threonine-protein kinase Nek5 protein (372-468AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase Nek5 (EC 2.7.11.1) (Never in mitosis A-related kinase 5) (NimA-related protein kinase 5), NEK5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NEK5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.