



G2E3 Antibody, FITC conjugated

Product Code	CSB-PA741972LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L622
Immunogen	Recombinant Human G2/M phase-specific E3 ubiquitin-protein ligase protein (240-337AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G2/M phase-specific E3 ubiquitin-protein ligase (EC 2.3.2.26) (G2/M phase-specific HECT-type E3 ubiquitin transferase), G2E3, KIAA1333
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	G2E3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.