



ENOSF1 Antibody, HRP conjugated

| | |
|---------------------|---|
| Product Code | CSB-PA748723LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q7L5Y1 |
| Immunogen | Recombinant Human Mitochondrial enolase superfamily member 1 protein (171-286AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Mitochondrial enolase superfamily member 1 (EC 4.2.1.68) (Antisense RNA to thymidylate synthase) (rTS) (L-fuconate dehydratase), ENOSF1, RTS TYMSAS |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | ENOSF1 |