



ENOSF1 Antibody, Biotin conjugated

Product Code	CSB-PA748723LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L5Y1
Immunogen	Recombinant Human Mitochondrial enolase superfamily member 1 protein (171-286AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial enolase superfamily member 1 (EC 4.2.1.68) (Antisense RNA to thymidylate synthase) (rTS) (L-fuconate dehydratase), ENOSF1, RTS TYMSAS
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ENOSF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.