



OTUD7B Antibody, HRP conjugated

Product Code	CSB-PA753588LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6GQQ9
Immunogen	Recombinant Human OTU domain-containing protein 7B protein (235-347AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	OTU domain-containing protein 7B (EC 3.4.19.12) (Cellular zinc finger anti-NF-kappa-B protein) (Cezanne) (Zinc finger A20 domain-containing protein 1) (Zinc finger protein Cezanne), OTUD7B, ZA20D1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	OTUD7B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.