



# FASTKD5 Antibody, HRP conjugated

Product Code	CSB-PA759173LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L8L6
Immunogen	Recombinant Human FAST kinase domain-containing protein 5, mitochondrial protein (516-664AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FAST kinase domain-containing protein 5, mitochondrial, FASTKD5, KIAA1792
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FASTKD5