



FASTKD5 Antibody, Biotin conjugated

Product Code	CSB-PA759173LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L8L6
Immunogen	Recombinant Human FAST kinase domain-containing protein 5, mitochondrial protein (516-664AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FAST kinase domain-containing protein 5, mitochondrial, FASTKD5, KIAA1792
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FASTKD5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.