



DODA Antibody

Product Code	CSB-PA767105LA01PIAD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7XA48
Immunogen	Recombinant Portulaca grandiflora 4,5-DOPA dioxygenase extradiol protein (1-271AA)
Raised In	Rabbit
Species Reactivity	Portulaca grandiflora
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	4,5-DOPA dioxygenase extradiol (EC 1.13.11.29), DODA
Immunogen Species	Portulaca grandiflora (Rose moss)
Research Area	Others
Target Names	DODA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.