



PEX26 Antibody, Biotin conjugated

Product Code	CSB-PA767196LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z412
Immunogen	Recombinant Human Peroxisome assembly protein 26 protein (1-159AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisome assembly protein 26 (Peroxin-26), PEX26
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PEX26
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.