



NLRC5 Antibody, FITC conjugated

Product Code	CSB-PA771457LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WI3
Immunogen	Recombinant Human Protein NLRC5 protein (95-208AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein NLRC5 (Caterpillar protein 16.1) (CLR16.1) (Nucleotide-binding oligomerization domain protein 27) (Nucleotide-binding oligomerization domain protein 4), NLRC5, NOD27 NOD4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NLRC5