



NLRC5 Antibody, Biotin conjugated

Product Code	CSB-PA771457LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WI3
Immunogen	Recombinant Human Protein NLRC5 protein (95-208AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein NLRC5 (Caterpillar protein 16.1) (CLR16.1) (Nucleotide-binding oligomerization domain protein 27) (Nucleotide-binding oligomerization domain protein 4), NLRC5, NOD27 NOD4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NLRC5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.