




PIGW Antibody

Product Code	CSB-PA773605LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z7B1
Immunogen	Recombinant Human Phosphatidylinositol-glycan biosynthesis class W protein (392-448AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol-glycan biosynthesis class W protein (PIG-W) (EC 2.3.-.-), PIGW
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PIGW

Image



Western Blot

Positive WB detected in: Rat kidney tissue

All lanes: PIGW antibody at 3.7μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 57 kDa

Observed band size: 57 kDa