



# PIGW Antibody, Biotin conjugated

|                     |  |
|---------------------|--|
| Product Code        | CSB-PA773605LD01HU   |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                          |
| Uniprot No.         | Q7Z7B1   |
| Immunogen           | Recombinant Human Phosphatidylinositol-glycan biosynthesis class W protein (392-448AA) |
| Raised In           | Rabbit   |
| Species Reactivity  | Human  |
| Tested Applications | ELISA  |
| Form                | Liquid   |
| Conjugate           | Biotin   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4       |
| Purification Method | >95%, Protein G purified   |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Alias               | Phosphatidylinositol-glycan biosynthesis class W protein (PIG-W) (EC 2.3.-.-), PIGW    |
| Immunogen Species   | Homo sapiens (Human)   |
| Research Area       | Tags & Cell Markers  |
| Target Names        | PIGW   |