



PNPLA5 Antibody, HRP conjugated

Product Code	CSB-PA773803LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z6Z6
Immunogen	Recombinant Human Patatin-like phospholipase domain-containing protein 5 protein (90-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Patatin-like phospholipase domain-containing protein 5 (EC 3.1.1.-) (GS2-like protein), PNPLA5, GS2L
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PNPLA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.