



LAPTM4B Antibody, FITC conjugated

Product Code	CSB-PA774814LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86VI4
Immunogen	Recombinant Human Lysosomal-associated transmembrane protein 4B protein (174-226AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysosomal-associated transmembrane protein 4B (Lysosome-associated transmembrane protein 4-beta), LAPTM4B
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	LAPTM4B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.