



CSMD3 Antibody, HRP conjugated

Product Code	CSB-PA800102LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z407
Immunogen	Recombinant Human CUB and sushi domain-containing protein 3 protein (233-380AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CUB and sushi domain-containing protein 3 (CUB and sushi multiple domains protein 3), CSMD3, KIAA1894
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CSMD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.