



# CSMD3 Antibody, Biotin conjugated

Product Code	CSB-PA800102LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z407
Immunogen	Recombinant Human CUB and sushi domain-containing protein 3 protein (233-380AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CUB and sushi domain-containing protein 3 (CUB and sushi multiple domains protein 3), CSMD3, KIAA1894
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CSMD3