



WDFY3 Antibody, FITC conjugated

Product Code	CSB-PA809030LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZQ1
Immunogen	Recombinant Human WD repeat and FYVE domain-containing protein 3 protein (1-100AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat and FYVE domain-containing protein 3 (Autophagy-linked FYVE protein) (Alfy), WDFY3, KIAA0993
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	WDFY3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.