




SMC5 Antibody

Product Code	CSB-PA810294LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IY18
Immunogen	Recombinant Human Structural maintenance of chromosomes protein 5 protein (296-403AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Structural maintenance of chromosomes protein 5 (SMC protein 5) (SMC-5) (hSMC5), SMC5, KIAA0594 SMC5L1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SMC5

Image



Western Blot

Positive WB detected in: HepG2 whole cell lysate, 293T whole cell lysate

All lanes: SMC5 antibody at 4.6µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 129 kDa

Observed band size: 145 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.