



SMC5 Antibody, Biotin conjugated

Product Code	CSB-PA810294LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IY18
Immunogen	Recombinant Human Structural maintenance of chromosomes protein 5 protein (296-403AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Structural maintenance of chromosomes protein 5 (SMC protein 5) (SMC-5) (hSMC5), SMC5, KIAA0594 SMC5L1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SMC5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.