



CERS3 Antibody, Biotin conjugated

Product Code	CSB-PA811597LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IU89
Immunogen	Recombinant Human Ceramide synthase 3 protein (59-134AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ceramide synthase 3 (CerS3) (EC 2.3.1.24) (Dihydroceramide synthase 3) (LAG1 longevity assurance homolog 3), CERS3, LASS3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CERS3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.