



FAM20C Antibody, Biotin conjugated

Product Code	CSB-PA816901LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXL6
Immunogen	Recombinant Human Extracellular serine/threonine protein kinase FAM20C protein (107-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Extracellular serine/threonine protein kinase FAM20C (EC 2.7.11.1) (Dentin matrix protein 4) (DMP-4) (Golgi casein kinase) (Golgi-enriched fraction casein kinase) (GEF-CK), FAM20C, DMP4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FAM20C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.