



FBXL20 Antibody, Biotin conjugated

Product Code	CSB-PA822228LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96IG2
Immunogen	Recombinant Human F-box/LRR-repeat protein 20 protein (1-137AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box/LRR-repeat protein 20 (F-box and leucine-rich repeat protein 20) (F-box/LRR-repeat protein 2-like), FBXL20, FBL2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FBXL20
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.