



PASK Antibody, HRP conjugated

Product Code	CSB-PA822296LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RG2
Immunogen	Recombinant Human PAS domain-containing serine/threonine-protein kinase protein (482-602AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PAS domain-containing serine/threonine-protein kinase (PAS-kinase) (PASKIN) (hPASK) (EC 2.7.11.1), PASK, KIAA0135
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PASK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.