



SAT2 Antibody, HRP conjugated

Product Code	CSB-PA836208LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96F10
Immunogen	Recombinant Human Diamine acetyltransferase 2 protein (1-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diamine acetyltransferase 2 (EC 2.3.1.57) (Polyamine N-acetyltransferase 2) (Spermidine/spermine N(1)-acetyltransferase 2) (Thialysine N-epsilon-acetyltransferase), SAT2, SSAT2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SAT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.