



NEK1 Antibody, Biotin conjugated

Product Code	CSB-PA836281LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PY6
Immunogen	Recombinant Human Serine/threonine-protein kinase Nek1 protein (626-809AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase Nek1 (EC 2.7.11.1) (Never in mitosis A-related kinase 1) (NimA-related protein kinase 1) (Renal carcinoma antigen NY-REN-55), NEK1, KIAA1901
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NEK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.