



SLC5A8 Antibody, Biotin conjugated

Product Code	CSB-PA836676LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N695
Immunogen	Recombinant Human Sodium-coupled monocarboxylate transporter 1 protein (540-610AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-coupled monocarboxylate transporter 1 (Apical iodide transporter) (Electrogenic sodium monocarboxylate cotransporter) (Sodium iodide-related cotransporter) (Solute carrier family 5 member 8), SLC5A8, AIT SMCT SMCT1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SLC5A8