



LILRA3 Antibody, Biotin conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA836677LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8N6C8 |
| Immunogen | Recombinant Human Leukocyte immunoglobulin-like receptor subfamily A member 3 protein (222-313AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Leukocyte immunoglobulin-like receptor subfamily A member 3 (CD85 antigen-like family member E) (Immunoglobulin-like transcript 6) (ILT-6) (Leukocyte immunoglobulin-like receptor 4) (LIR-4) (Monocyte inhibitory receptor HM43/HM31) (CD antigen CD85e), LILRA3, ILT6 LIR4 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | LILRA3 |