



ACD Antibody, FITC conjugated

Product Code	CSB-PA839289LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AP0
Immunogen	Recombinant Human Adrenocortical dysplasia protein homolog protein (328-497AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adrenocortical dysplasia protein homolog (POT1 and TIN2-interacting protein), ACD, PIP1 PTOP TINT1 TPP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ACD