



BMPER Antibody, Biotin conjugated

Product Code	CSB-PA839829LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N8U9
Immunogen	Recombinant Human BMP-binding endothelial regulator protein (271-418AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BMP-binding endothelial regulator protein (Bone morphogenetic protein-binding endothelial cell precursor-derived regulator) (Protein crossveinless-2) (hCV2), BMPER, KIAA1965
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BMPER
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.