



NETO1 Antibody, Biotin conjugated

Product Code	CSB-PA840571LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDF5
Immunogen	Recombinant Human Neuropilin and tolloid-like protein 1 protein (367-507AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuropilin and tolloid-like protein 1 (Brain-specific transmembrane protein containing 2 CUB and 1 LDL-receptor class A domains protein 1), NETO1, BTCL1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NETO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.