



ELMOD3 Antibody, FITC conjugated

Product Code	CSB-PA842670LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96FG2
Immunogen	Recombinant Human ELMO domain-containing protein 3 protein (1-102AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ELMO domain-containing protein 3 (RNA-binding motif and ELMO domain-containing protein 1) (RNA-binding motif protein 29) (RNA-binding protein 29), ELMOD3, RBED1 RBM29
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ELMOD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.