



SIGLEC10 Antibody, Biotin conjugated

Product Code	CSB-PA842706LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LC7
Immunogen	Recombinant Human Sialic acid-binding Ig-like lectin 10 protein (572-697AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sialic acid-binding Ig-like lectin 10 (Siglec-10) (Siglec-like protein 2), SIGLEC10, SLG2
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SIGLEC10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.