



FAM84A Antibody, HRP conjugated

Product Code	CSB-PA846652LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96KN4
Immunogen	Recombinant Human Protein FAM84A protein (15-144AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM84A (Neurologic sensory protein 1) (NSE1), FAM84A, NSE1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM84A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.