



KCNV2 Antibody, Biotin conjugated

Product Code	CSB-PA848410LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDN2
Immunogen	Recombinant Human Potassium voltage-gated channel subfamily V member 2 protein (1-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium voltage-gated channel subfamily V member 2 (Voltage-gated potassium channel subunit Kv8.2), KCNV2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNV2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.