



LIPH Antibody, FITC conjugated

Product Code	CSB-PA848833LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWY8
Immunogen	Recombinant Human Lipase member H protein (175-325AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lipase member H (LIPH) (EC 3.1.1.-) (LPD lipase-related protein) (Membrane-associated phosphatidic acid-selective phospholipase A1-alpha) (mPA-PLA1 alpha) (Phospholipase A1 member B), LIPH, LPDLR MPAPLA1 PLA1B
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LIPH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.