



PIK3IP1 Antibody, Biotin conjugated

Product Code	CSB-PA850293LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96FE7
Immunogen	Recombinant Human Phosphoinositide-3-kinase-interacting protein 1 protein (22-168AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphoinositide-3-kinase-interacting protein 1 (Kringle domain-containing protein HGFL), PIK3IP1, HGFL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PIK3IP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.