



CDCA7L Antibody, FITC conjugated

Product Code	CSB-PA850301LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GN5
Immunogen	Recombinant Human Cell division cycle-associated 7-like protein (14-167AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell division cycle-associated 7-like protein (Protein JPO2) (Transcription factor RAM2), CDCA7L, HR1 JPO2 R1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CDCA7L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.