



CDCA7L Antibody, FITC conjugated

Product Code	CSB-PA850301LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GN5
Immunogen	Recombinant Human Cell division cycle-associated 7-like protein (14-167AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell division cycle-associated 7-like protein (Protein JPO2) (Transcription factor RAM2), CDCA7L, HR1 JPO2 R1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CDCA7L