



VPS52 Antibody, HRP conjugated

Product Code	CSB-PA850779LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N1B4
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 52 homolog protein (277-377AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 52 homolog (SAC2 suppressor of actin mutations 2-like protein), VPS52, SACM2L
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VPS52
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.