



NKD1 Antibody, Biotin conjugated

Product Code	CSB-PA853369LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969G9
Immunogen	Recombinant Human Protein naked cuticle homolog 1 protein (269-384AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein naked cuticle homolog 1 (Naked-1) (hNkd) (hNkd1), NKD1, NKD
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	NKD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.