



FERMT2 Antibody, Biotin conjugated

Product Code	CSB-PA853382LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AC1
Immunogen	Recombinant Human Fermitin family homolog 2 protein (155-232AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fermitin family homolog 2 (Kindlin-2) (Mitogen-inducible gene 2 protein) (MIG-2) (Pleckstrin homology domain-containing family C member 1) (PH domain-containing family C member 1), FERMT2, KIND2 MIG2 PLEKHC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FERMT2